


DMG MICROLABS


ATTENTION

Please read instructions before using Q-Swabs


DIRECTIONS FOR USE:



1. LABEL SWAB.



2. HOLD SWAB BY LABEL. TWIST TOP NEAR RED AREA AND REMOVE SWAB. DO NOT BREAK BULB CONTAINING FLUID.



3. USE DRY SWAB TO COVER A 100cm² AREA



4. REPLACE SWAB INTO CASE. DO NOT BREAK BULB CONTAINING FLUID.

5. COURIER TO LAB WITH ICEPACKS.



Accreditation
No.: 14944

1/21 Smallwood Place, Murarrie. QLD. 4172
Phone: 07 3902 0575
Fax: 07 3395 7275
ABN: 19 001 207 166

PO Box 3410 TINGALPA DC
Qld 4173 Australia
Email: admin@dmgmicrolabs.com.au
Web : www.dmgmicrolabs.com.au

Revised 24 July 2008